

Versa



An exam table as diverse as your
patient population.



*Standard top with patient assist
handles, stirrups and pelvic tilt*



Standard top



Pediatric top

Brewer

Versa

gives you the ability to readily and cost-effectively adapt to your changing needs. When it comes to care, the standard model provides a new standard of quality and efficiency for most practices. Easily convert to a pediatric top to safely and effectively care for children from newborns to teenagers.



IMPROVED PATIENT EXPERIENCE

- Stable cushioned platform with 500 pound patient weight capacity comfortably supports patients of all builds.
- Large, sturdy step with anti-slip surface enhances patient security, and is easy to clean.
- Multi-position patient assist handles offer security and support patient self-positioning.



HIGHLY CUSTOMIZABLE

- Easy to connect, removable upholstery top is available in standard adult and pediatric sizes to meet the changing needs of your facility.
- Versa is truly customizable – choose only the options you need to stay within your specific budget.
- Optional features include patient assist handles, electrical receptacle, drawer warmer and pelvic tilt, debris pan, urology drain pan, exam light bracket, knee crutch attachments and front drawer dividers.



MORE CONVENIENT ACCESSIBILITY

- Ergonomically designed pass-through side drawer positioned purposefully at waist height to reduce bending and twisting for caregivers while offering a generous 3.7 ft³ of storage space.
- Front drawer strategically positioned at lower leg height to eliminate uncomfortable situations for caregivers and patients.
- Rounded corners and softened edges throughout safeguard patients and caregivers while improving patient access.



Large pass-through drawer with included dividers offers 3.7 ft³ of storage to organize essentials.



Standard top

- Easy connect upholstery top can be removed for a deep cleaning or interchanged with a pediatric top.

- 500 pound patient weight capacity

- Ergonomically located side drawer reduces caregiver bending and twisting

Pediatric top

- Sloped and cushioned baby-cradle edges keep younger children safely in place during exams

Total Storage Capacity

Brewer Versa 1 Drawer Front; 1 Side	4.7 ft ³
Brewer Access 5000 2 Drawer Front; 2 Side	5.2 ft ³
Competitive Table 2 Drawer Front; 3 Side	4.2 ft ³

- Lower drawer positioned so caregivers do not have to infringe on patient personal space

- Large, anti-slip patient step

- Padded leg extension accommodates older children

Multi-position patient assist handles designed to expedite patient transfers and facilitate efficient repositioning.

Women's health package includes drawer warmer, pelvic tilt and ergonomic stirrups.

Large, sturdy step with integral anti-slip surface enhances patient safety and is easy to clean.

Versa

The added versatility you need to enhance care excellence across your patient population.

Options



Patient Assist Handles



Ergonomic Stirrups



Paper Fastener Strap



Receptacle & Heater

Accessories



Knee Crutches



Stainless Steel Debris Pan



Urology Drain Pan



Welch Allyn Bracket

Configure your exam table

Here's how to build your Versa exam table.

Top	Assist Handles	Stirrups	Paper Roll Holder	Pelvic Tilt	Receptacle & Heater	Power Cord	Upholstery Color
1000 <input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
AX- Standard PX- Pediatric	H- Yes X- No	S- Yes X- No	D- 18" or 21" rolls	T- Yes X- No	RH- Outlet & Heater RX- Outlet Only XX- Neither	X	*21 to 32 - Std Colors SP - Special Color

* Assist handles, electric & power cord are not available with pediatric top.

Example: 1000AX-HSDX-XXX-21 creates an order for a Versa with standard top, assist handles, stirrups and an 18 to 21-in. paper roll holder in clamshell upholstery.

Recommended Configurations

Base Model: 1000AX-XXDX-XXX- [2 digit color code]
 Pediatric Base Model: 1000PX-XXDX-XXX- [2 digit color code]
 Women's Health Base Model: 1000AX-XSDT-RHX- [2 digit color code]
 Standard model with assist handles: 1000AX-HXDX-XXX- [2 digit color code]

Available in 12 Standard Color Choices

Recommended _____ Also Available _____



24-Black Satin 29-Deep Sea 30-Ivy 31-Blue Fog 32-Deep Sapphire 21-Clamshell 22-Gunmetal 23-Feather 25-Saddle 26-Tapestry 27-Cabernet 28-Cocoa

Brewer

The Brewer Company, LLC
 N88 W13901 Main Street, Suite 100
 Menomonee Falls, WI 53051

1.888.273.9371



To discover all the ways you'll benefit from Versa's versatility, visit www.brewercompany.com

We Help Healers Heal.

LIT-M-Q55-REV C

Shop tables & accessories at