

# Brewer

## Total Storage Capacity

### Versa

1 Drawer Front; 1 Side 4.7 ft<sup>3</sup>

### Access 5000

2 Drawer Front; 2 Side 5.2 ft<sup>3</sup>

### Competitive Table

2 Drawer Front; 3 Side 4.2 ft<sup>3</sup>

\* Measurements are approximate



## Available in 12 standard colors

### Recommended



24-Black  
Satin

29-Deep Sea

30-Ivy

31-Blue Fog

32-Sapphire

### Also available



21-Clamshell

22-Gunmetal

23-Feather

25-Saddle

26-Tapestry

27-Cabernet

28-Cocoa

# Versa

## Exam Table



Leading-edge design & added versatility help you enhance care across patient populations

The Versa Exam Table readily and cost-effectively adapts to your changing patient and practice needs.



### HIGHLY CUSTOMIZABLE

- Easy-connect, removable upholstery top is available in standard & pediatric configurations
- Add options as needed to economically meet practice requirements



### IMPROVED PATIENT EXPERIENCE

- Large, sturdy step with integral, anti-slip surface enhances patient safety
- Stable, cushioned platform with 500-lb patient weight capacity comfortably supports patients of all builds
- Optional multi-position patient-assist handles sized for more confident transitions and easier patient self-positioning
- Seamless upholstery & integral step surface offer easy cleanability and infection control



### MORE ERGONOMIC ACCESS

- Rounded corners and softened edges throughout safeguard patients and caregivers, while improving patient access
- Pass-through side drawer purposefully positioned at waist height to reduce bending and twisting
- Front drawer strategically positioned so caregivers do not infringe on a patient's personal space

**We Help Healers Heal**

Contact Brewer at 1.888.273.9371 to get a quote or schedule an on-site demo.

## Product Specifications

- Weight capacity: 500 lbs. (226 kg)
- Seat height: 32.5" (82.55 cm)
- Length: 58"; extendable to 72" (147.3 cm-182.9 cm)
- Pediatric table length: 51"; expandable to 65" (129.5 cm-165.1 cm)
- Pediatric table width: 26" (66 cm)
- Leg extension length: 14" (35.5 cm)
- Upholstery width: 27.25" (69.2 cm)
- Front storage drawer: 21.3" x 19.75" x 4.25" (59.9 cm x 50.1 cm x 9.7 cm)
- Side storage drawer: 23.6" x 35.75" x 7.6" (57.5 cm x 90.8 cm x 19.3 cm)
- Paper roll holder: accommodates up to 21" paper rolls
- UL certified: IEC60601-1 Ed. 3.1

## Standard Features

- Industry-best 5-year warranty
- Easy-connect, removable upholstery top
- Large, sturdy patient step with integral anti-slip surface
- Easy-to-clean, seamless upholstery available in 12 standard color choices
- Numerous additional custom colors available to choose from
- Two side drawer dividers for organization

## Optional Features

- Hospital-grade electric receptacle
- Paper fastener strap on standard top offers convenience
- Patient assist handles
- 120VAC, 60Hz
- 60-Watt front drawer heater
- 5° pelvic tilt
- Ergonomic Stirrups
- Footrest pad

## Models

- Standard model:  
1000AX-XXDX-XXX-\_\_ [2 digit color code]
- Standard model with pelvic tilt, front drawer heater and outlet:  
1000AX-XSDT-RHX-\_\_ [2 digit color code]
- Standard model with assist handles:  
1000AX-HXDX-XXX-\_\_ [2 digit color code]
- Pediatric model:  
1000PX-XXDX-XXX-\_\_ [2 digit color code]

## Optional Accessories

- 98464 – Stainless Steel Treatment Pan
- 3500286 – Deluxe articulating knee crutches with assembly kit
- 100253 – Urology Drain Pan System
- 3500288 – Front Drawer Dividers
- 3500347 – Brewer Welch Allyn® Exam Light / IV Bracket (factory installed)
- 3500398 – Brewer Welch Allyn® Exam Light / IV Bracket (field installed)
- 3500237 – Paper Fastener Strap
- 3500077 – Footrest Pad

\* Product specification measurements are approximate and may vary +/- 0.5"



# Brewer

## We Help Healers Heal

Learn more at [brewercompany.com](http://brewercompany.com)

Shop tables & accessories at



LIT-M-056