

Flat Top, Classic Series, Straight Line Treatment Table



MODEL 100

- Solid hardwood legs
- H-brace construction
- Legs have pre-inserted carriage bolts for fast easy assembly
- Seamless, rounded corner top
- 2" of firm foam padding (5 cm)
- Available in Natural and Dark Cherry finish
- Choice of upholstery colors
- Fast, easy assembly with basic tools
- 400 lbs. load capacity under normal use (181.43 kg)

SPECIFICATIONS

| | |
|--------|-----------------------|
| 100-24 | 72" L x 24" W x 31" H |
| 100-27 | 72" L x 27" W x 31" H |
| 100-30 | 72" L x 30" W x 31" H |
| 100-38 | 78" L x 30" W x 31" H |

California residents please be advised of the following, pursuant to Proposition 65:

WARNING: This product can expose you to chemicals, including diisononyl phthalate (DINP), which is known to the State of California to cause cancer. For more information go to www.P65Warnings.ca.gov.

2

Options



Finishes



Upholstery Colors

