

## Classic Series, Straight Line Treatment Table with Drawers



### MODEL 1013

- Features hardwood legs, laminate shelf, two easy clean, plastic drawers and adjustable backrest
- Backrest has 6 settings of approximately: 60°, 55°, 50°, 43°, 36° & 27°
- Drawers can be positioned at either end, and from either side
- Scratch resistant topcoat finish
- Fast and easy assembly
- Seamless rounded corner top
- 2" firm foam padding (5 cm)
- Heavyweight, knit-backed upholstery
- 400 lbs. load capacity under normal use (181.43 kg)

### SPECIFICATIONS

1013-27	72" L x 27" W x 31" H
1013-30	72" L x 30" W x 31" H
1013-38	78" L x 30" W x 31" H

### Prop 65 Warning

**WARNING:** This product can expose you to chemicals, including diisononyl phthalate (DINP), which is known to the State of California to cause cancer. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

2

### Options



20- small Pillow

077- Positioning Strap

034- Paper Dispenser & Cutter Combo

040- Paper Cutter

030- Paper Dispenser

005- Special Height

006- Special Length


009- Nose Hole

055- Drawer Locks

### Finishes



Natural



Dark Cherry

### Upholstery Colors



Slate Blue



Wedgewood



Royal Blue



Desert Tan



Burgundy



Country Mist



Warm Gray



Gunmetal



Black