

## Flat Top, Classic Series, Straight Line Treatment Table with Shelf



### MODEL 200

- Solid hardwood legs
- H-brace construction
- Full laminate shelf in Natural or Dark Cherry to coordinate with frame
- Legs have pre-inserted carriage bolts for fast easy assembly
- Seamless, rounded corner top
- 2" of firm foam padding (5 cm)
- Available in Natural and Dark Cherry finish
- Choice of upholstery colors
- Fast, easy assembly with basic tools
- 400 lbs. load capacity under normal use (181.43 kg)

### SPECIFICATIONS

200-24	72" L x 24" W x 31" H
200-27	72" L x 27" W x 31" H
200-30	72" L x 30" W x 31" H
200-38	78" L x 30" W x 31" H


California residents please be advised of the following, pursuant to Proposition 65:



**WARNING:** This product can expose you to chemicals, including diisononyl phthalate (DINP), which is known to the State of California to cause cancer. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

2

### Options



20- small Pillow



077- Positioning Strap



034- Paper Dispenser & Cutter Combo



040- Paper Cutter




030- Paper Dispenser



005-Special Height



006-Special Length



009-Nose Hole

### Finishes



Natural



Dark Cherry

### Upholstery Colors



Slate Blue



Wedgewood



Royal Blue



Desert Tan



Burgundy



Country Mist



Warm Gray



Gunmetal



Black