

## Alpha Series, Treatment Table with Shelf




### MODEL 3020

- Features powder-coated, all-steel frame, gray laminate shelf and an adjustable backrest
- Easy bolt-through assembly
- Paper dispenser included
- Adjustable leg levelers
- Seamless rounded corner top
- 2" firm foam padding (5 cm)
- Heavyweight, knit-backed upholstery
- 450 lbs. load capacity under normal use (204 kg)


### SPECIFICATIONS

3020-27	72" L x 27" W x 31" H
302030	72" L x 30" W x 31" H

California residents please be advised of the following, pursuant to Proposition 65:

 **WARNING:** This product can expose you to chemicals, including diisononyl phthalate (DINP), which is known to the State of California to cause cancer. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov). **2**

### Options



### Upholstery Colors

