

## Alpha Series, Treatment Table with Shelving



### MODEL 3030

- Features powder-coated, all-steel frame, gray laminate shelving and an adjustable backrest
- Easy bolt-through assembly
- Paper dispenser included
- Adjustable leg levelers
- Seamless rounded corner top
- 2" firm foam padding (5 cm)
- Heavyweight, knit-backed upholstery
- 450 lbs. load capacity under normal use (204 kg)

### SPECIFICATIONS


3030-27	72" L x 27" W x 31" H
3030-30	72" L x 30" W x 31" H

California residents please be advised of the following, pursuant to Proposition 65:

**!** **WARNING:** This product can expose you to chemicals, including diisononyl phthalate (DINP), which is known to the State of California to cause cancer. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

2

### Options



020 Small Pillow

060 Large Pillow

009 Nosehole

077 Safety Strap

040 Paper Cutter

### Upholstery Colors



Slate Blue

Wedgewood

Royal Blue

Desert Tan

Burgundy

Country Mist

Warm Gray

Gunmetal

Black