

## Flat Top, Alpha Series, Treatment Table with Drawers



### MODEL 3130

- Gray, powder-coated, all-welded, steel frame
- Gray laminate drawer fronts and shelf
- Adjustable height levelers
- Pre-threaded hardware for easy assembly
- Paper dispenser included
- 2" of firm foam padding (5 cm)
- Traditional square corner top
- Fast, easy assembly with basic tools
- 450 lbs. load capacity under normal use (204 kg)

### SPECIFICATIONS

3130-27	72" L x 27" W x 31" H
3130-30	72" L x 30" W x 31" H

California residents please be advised of the following, pursuant to Proposition 65:

**!** **WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

1

### Options



020 Small Pillow

060 Large Pillow

009 Nosehole

077 Safety Strap

040 Paper Cutter

055 Drawer Locks

### Upholstery Colors



Slate Blue

Wedgewood

Royal Blue

Desert Tan

Burgundy

Country Mist

Warm Gray

Gunmetal

Black