

## Lab X Series, Extra-Wide, Blood Drawing Chair with Padded Arms

### MODEL 66000

- Padded seat and backrest
- Height and depth adjustable, dual, padded, flip arms
- Flip arm(s) positions in place when in use, then flips out of the way for easy exit
- Dual ClintonClean elbow rests
- ClintonClean™ elbow rests
- A full 35" wide seat to accommodate a larger patient, or parent and child (89 cm)
- 1-1/4" square, heavy duty, all-welded, tubular steel frame (3.18 cm)
- Comfortable 2" padding throughout (5 cm)
- Neutral Gray, powder-coated frame
- Adjustable foot levelers
- Easy-clean, durable vinyl in a selection of colors
- ClintonClean™ arms option available
- 400 lbs. combined load capacity under normal use (181.4 kg)



**IMPORTANT NOTE:** Clinton Blood Drawing Chairs are not engineered to restrain a fainting or falling occupant. The chair's occupant should NEVER be left unattended.

**ORDERING NOTE:** MUST add arm suffix number to model number when ordering any model with arm option, otherwise chair will be shipped as shown.

**All blood drawing products are non-returnable unless they are unused in unopened, original packaging.**



California residents please be advised of the following, pursuant to Proposition 65:

**WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

### MODEL 66000 SPECIFICATIONS

Seat	35" W x 18" D x 20" H
Overall	46" W x 26" D x 36" H
Arm Height Range	27" — 33"
Arm Adjustment	6"

### OPTIONS



665 ClintonClean™, Armrests and Flip Arm(s)

### Standard Upholstery Colors

