

Lab X Series, Bariatric, Blood Drawing Chair with Padded Arms

MODEL 66000B

- Padded, true, bariatric seat features extra width and depth for added comfort
- Comfortable padded arched backrest
- Height and depth adjustable, dual, padded, flip arms
- Flip arm(s) positions in place when in use, then flips out of the way for easy exit
- Dual ClintonClean™ arm rests
- A full 35" wide seat (89 cm)
- 1-1/4" square, heavy duty, all-welded, tubular steel frame (3.18 cm)
- Comfortable 2" padding throughout (5 cm)
- Neutral Gray powder-coated frame
- Easy-clean, durable vinyl in a selection of colors
- ClintonClean™ arms option available
- 700 lbs. combined load capacity under normal use (317.8 kg)



IMPORTANT NOTE: Clinton Blood Drawing Chairs are not engineered to restrain a fainting or falling occupant. The chair's occupant should NEVER be left unattended.

ORDERING NOTE: MUST add arm suffix number to model number when ordering any model with arm option, otherwise chair will be shipped as shown.

All blood drawing products are non-returnable unless they are unused in unopened, original packaging.

MODEL 66000B SPECIFICATIONS

Seat	35" W x 21" D x 20" H
Overall	46" W x 29" D x 36" H
Arm Height Range	27" — 33"
Arm Adjustment	6"




OPTIONS



665 ClintonClean™, Armrests and Flip Arm(s)

California residents please be advised of the following, pursuant to Proposition 65:

1
 **WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

Standard Upholstery Colors

