

Select Series, Pediatric Scale Table

MODEL 7830

TABLE FEATURES

- Comfortable work height of 36"
- Top is rimmed on three sides for safety
- Top's front and back rim edges are beveled and the corners rounded for safety
- Easy-clean, replaceable center pad
- Built-in infantometer with folding ends for measuring height
- Easy access to weight, measurement and storage
- One drawer and two-door storage
- Built-in concealed paper dispenser and cutter
- Recessed front toe kick
- Easy-clean all laminate base

SCALE FEATURES

- Dependable, accurate Doran® scale
- Easy-read digital display
- Auto calibration on start-up
- Auto shut-off
- Electric or battery operation (6 AA batteries required)
- Weights in increments of lbs., ozs. or kilos
- Minimum registered weight is 5 lbs.—Maximum weight for scale is 40 lbs.



MODEL 7830 SPECIFICATIONS

43"W x 24"D x 38.5"H

California residents please be advised of the following, pursuant to Proposition 65:

WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

1



Options



066 Door & Inside
Latch Combo



055 Drawer Locks



078-S
Seal Scale Pal



078-B
Bear Scale Pal



078-K
Kitty Scale Pal



078-D
Doggie Scale Pal




Select Series
Graphics

Laminate Colors



Maple



Dark Cherry



Gray



Slate Gray

Standard Upholstery



Slate Blue



Wedgewood



Royal Blue



Desert Tan



Burgundy



Country Mist



Warm Gray



Gunmetal



Black

Features




Infantometer

Built-in infantometer with folding ends for measuring height



Paper Dispenser

Choose to access concealed paper dispenser from either end



Rimmed Top with Beveled Edges

Top is rimmed on three sides for safety and the front and back rim edges are beveled and the corners rounded for safety



Padded Top

Easy-clean, padded, upholstered top for comfort



Scale Head

Easy-read, digital display