

## General Ultrasound Table with Adjustable Backrest



### MODEL 80062

- Two-piece top with easy-to-adjust, pneumatic backrest
- Trendelenburg/Reverse Trendelenburg plus or minus 15 °
- Easy-to-use hand control
- All-welded base and top frame
- Unobstructed base design
- Dual column lift system adds rigidity
- Paper dispenser included
- Industry leading 5-year, limited warranty
- Optional side rails & casters available
- 34" wide model available
- Model 80062: 500 lbs load capacity under normal use
- Model 80062-X: 600 lbs load capacity under normal use



*Trendelenburg*



*Reverse Trendelenburg*



*Optional Side Rails & Casters*

### 80062 SPECIFICATIONS

	length	width	height range
overall	72"	27"	21 <sup>1</sup> / <sub>2</sub> " – 37 <sup>1</sup> / <sub>2</sub> "
backrest	38"	27"	

**80062-X** Extra-Wide 34" width model

California residents please be advised of the following, pursuant to Proposition 65:

**WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

1



### Upholstery Colors



*Slate Blue*

*Wedgewood*

*Royal Blue*

*Desert Tan*

*Burgundy*

*Country Mist*

*Warm Gray*

*Gunmetal*

*Black*

## OPTIONS

### Positioning Straps Set

- 2" x 92" strap is secured with quick hook and loop Velcro® fasteners
- 2 straps per set

080	length	width
	92"	2"
	233.67 cm	5 cm



### Safety Rails (pair)

- 3 height positions
- Plunger height adjustment

098



### Locking Casters

- Heavy duty, rubber wheel, easy-roll casters

087



### Paper Cutter

- Clear plastic with elastic strap and buckle

040

