

## Shrouded, PowerXL Table with Adjustable Backrest and Drop Section

### MODEL 81399

- Features easy-clean, 30" wide top with adjustable backrest and drop section
- 600 lbs. (272 kg) load capacity under normal use
- Easy-clean, ABS plastic enclosed lift
- Backrest adjusts from either side and raises 80 degrees from flat
- Drop section adjusts from either side and lowers up to 80 degrees from flat
- Smooth, electric, raise and lower height adjustment
- Pneumatic backrest cylinder is mounted internally
- Full width backrest hinge adds rigidity for greater support
- Full width backrest hinge adds rigidity for greater support
- Heavy duty, steel frame is fabricated in the USA by Clinton Industries
- Dependable scissors-lift features an H-leg design for extra stability
- Rounded corner 3/4" plywood top with 2" of ultra firm padding covered by heavy weight upholstery
- Harden steel bolts and self lubricating pivot

CLINTON OPEN BASE & SHROUDED  
**POWERXL**  
*Extra width, Extra capacity, Extra Affordable*



### 81399 SPECIFICATIONS

Overall	76" L x 30" W x 21-34" H
Backrest	38" L x 30" W
Drop Section	13" L x 30" W

bearings for years of dependable service

- Low height may help achieve ADA tax credit



**California residents please be advised of the following, pursuant to Proposition 65:**

**!** **WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

### Upholstery Colors



## OPTIONS

### Positioning Straps Set

- 2" x 92" strap is secured with quick hook and loop Velcro® fasteners
- 2 straps per set (92" x 2")

080



### Safety Rails (pair)

- 3 height positions
- Plunger height adjustment

082



### Hand-Held Controller

- Alternative to foot control
- Stores easily under top when not in use

084



### 030 Paper Dispenser


- Accommodates paper rolls from 18" (45.72cm) to 21" (53.34cm) wide

030

### Paper Cutter

- Clear plastic with elastic strap and buckle

040



### Paper Dispenser/Cutter Combo

- Paper dispenser and paper holder ordered as one option

034