

## Open Base, Extra Wide, Bariatric, Straight Top Power Table

### MODEL 84100

- Extra wide, bariatric model with choice of two widths
- Features easy-clean, one-piece top
- Electric power height adjustment with foot control
- Ultra-strong, uni-frame base and top construction
- Heavy duty, all welded, lift mechanism
- Durable gray powder-coated steel frame
- Self-lubricating pivot points
- Smooth 110 volt electric height adjustment
- Power supply and motor meets UL 60601-1 and CAN/CSA C22.2 601.1
- Hands-free low voltage foot control
- Easy service, self contained, electrical components
- Rounded corner, upholstered top
- 2" ( 5 cm) ultra-firm padding
- Low height may help achieve ADA tax credit
- Customize with options
- Made in the USA by Clinton Industries
- 600 lbs. (272 kg) load capacity under normal use

# POWER 600 BARIATRIC



### 84100-34 SPECIFICATIONS

Overall 72" L x 34" W x 18-35" H

84100-40 Extra-Wide 40" width model

California residents please be advised of the following, pursuant to Proposition 65:

**WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

1



### Upholstery Colors



Slate Blue



Wedgewood



Royal Blue



Desert Tan



Burgundy



Country Mist



Warm Gray



Gunmetal



Black

## OPTIONS

### Positioning Straps Set

- 2" x 92" strap is secured with quick hook and loop Velcro® fasteners
- 2 straps per set (92" x 2")

080

080



### Safety Rails (pair)

- 3 height positions
- Plunger height adjustment

082

082



### Hand-Held Controller

- Alternative to foot control
- Stores easily under top when not in use

084

084



### Locking Casters

- Heavy duty, rubber wheel, easy-roll casters

087

087



### 030 Paper Dispenser

- Accommodates paper rolls from 18" (45.72cm) to 21" (53.34cm) wide

030

030




### Paper Cutter

- Clear plastic with elastic strap and buckle

040

040



### Paper Dispenser/Cutter Combo

- Paper dispenser and paper cutter ordered as one option

034