

Open Base, Multi-Use Power Imaging Table with Stirrups

MODEL 85309

- Pneumatic backrest is adjustable from either side and raises 80 degrees from flat
- Footrest is adjustable from either side and lowers up to 80 degrees from flat
- 10" x 10" hinged imaging window section with metal slide latch
- Retractable, 2-position, length adjustable stirrups
- Electric power height adjustment with foot control
- Ultra-strong, uni-frame base and top construction
- Heavy duty, all welded, lift mechanism
- Durable gray powder-coated steel frame
- Self-lubricating pivot points
- Smooth 110 volt electric height adjustment
- Power supply and motor meets UL 60601-1 and CAN/CSA C22.2 601.1
- Hands-free low voltage foot control
- Easy service, self contained, electrical components
- Rounded corner, upholstered top
- 2" ultra-firm padding
- Low height may help achieve ADA tax credit
- Customize with options
- Made in the USA by Clinton Industries



MODEL 85309 SPECIFICATIONS

Overall	76" L x 27" W x 18-35" H
Backrest	38" L x 27" W
Footrest	13" L x 19" W
Drop Window	10" L x 10" W

California residents please be advised of the following, pursuant to Proposition 65:

1
WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.



Upholstery Colors



OPTIONS

Positioning Straps Set

- 2" x 92" strap is secured with quick hook and loop Velcro® fasteners
- 2 straps per set (92" x 2")

080

080



Safety Rails (pair)

- 3 height positions
- Plunger height adjustment

082

082



Hand-Held Controller

- Alternative to foot control
- Stores easily under top when not in use

084

084



Locking Casters

- Heavy duty, rubber wheel, easy-roll casters

087

087



030 Paper Dispenser

- Accommodates paper rolls from 18" (45.72cm) to 21" (53.34cm) wide

030

030




Paper Cutter

- Clear plastic with elastic strap and buckle

040

040



Paper Dispenser/Cutter Combo

- Paper dispenser and paper cutter ordered as one option

034