

## E-Series, Space Saver, Power Table with Drop Section

### MODEL 86220

- Drop Section is adjustable from either side and raises 80 degrees from flat.
- Standard drop section, but it may be ordered with a nose hole as an option (Please specify when ordering)
- Electric power height adjustment with foot control
- Ultra-strong, uni-frame base and top construction
- Heavy duty, all welded, lift mechanism
- Durable gray powder-coated steel frame
- Self-lubricating pivot points
- Smooth 110 volt electric height adjustment
- Power supply and motor meets UL 60601-1 and CAN/CSA C22.2 601.1
- Hands-free low voltage foot control
- Easy service, self contained, electrical components
- Seamless, rounded corner, upholstered top
- 1.5" (4 cm) ultra-firm padding
- Low height may help achieve ADA tax credit
- Customize with options
- Made in the USA by Clinton Industries
- 450 lbs. (204 kg) load capacity under normal use

**E-Series**  
POWER TABLES



#### SPECIFICATIONS

Overall 76" L x 27" W x 18.5 –35.5" H

Drop Section 13" L x 27" W

**California residents please be advised of the following, pursuant to Proposition 65:**

**!** **WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).



### Options



**009 Nose Hole** — Nose hole in drop section of 86220 or 86300

**030 Paper Dispenser** — Accommodates paper rolls from 18" (45.72cm) to 21" (53.34cm) wide

**040 Paper Cutter** — Clear plastic with elastic strap and buckle

**034 Paper Dispenser/Holder Combo** — Paper dispenser and paper holder ordered as one option

### Upholstery Colors

