

Cabinet Style, Space Saver Table with Step Stool

MODEL 8844

- Premium, vinyl upholstery over 2" firm, foam padding (5 cm)
- Full 27" width for patient comfort (69 cm)
- Pneumatic backrest adjusts from either side
- Pull out step stool with safety tread
- Pull-out, slate gray, laminate leg rest on steel track
- Easy care, Gray, laminate base
- Color coordinated doors and top available
- Chip resistant edges
- 1 large side storage compartment with shelf
- Storage compartment is accessible from both sides (no right or left handed tables)
- 2 flat panel drawers with all steel slides
- Built with plywood and laminate panels that are California CARB II compliant
- Top sections are easily replaceable in the field
- Adjustable feet for leveling
- Paper dispenser included
- 400 lbs. load capacity under normal use (181 kg)



8844 SPECIFICATIONS

8844 58"/72" extended L x 28.5" W x 31" H




California residents please be advised of the following, pursuant to Proposition 65:

WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.


1

Options



040- Paper Holder

Laminate Color



Gray

Upholstery Colors



Slate Blue

Wedgewood

Royal Blue

Desert Tan

Burgundy

Country Mist

Warm Gray

Gunmetal

Black