

Family Practice Table

MODEL 8870

- Premium, vinyl upholstery over 2" firm, foam padding (5 cm)
- Full 28.5" width for patient comfort (72.39 cm)
- Air-spring adjustable backrest
- Coated, two position adjustable stirrups
- Pull-out leg rest on steel track
- Easy care, Gray, laminate base
- Color coordinated doors and top available
- Chip resistant edges
- 1 large side storage compartment with shelf
- Storage compartment is accessible from both sides (no right or left handed tables)
- 2 drawers with all steel slides
- Built with plywood and laminate panels that are California CARB II compliant
- Top sections are easily replaceable in the field
- Adjustable feet for leveling
- Paper dispenser included
- 400 lbs. load capacity under normal use (181.43 kg)



SPECIFICATIONS

58"/72" extended L x 28.5" W x 31" H



California residents please be advised of the following, pursuant to Proposition 65:

! **WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

Features




Pull-out, slate, laminate, leg rest



2 drawer storage



2-Position Stirrups



Paper Roll Storage



Airlift Backrest



Both Sides Access

Laminate Color



Gray

Option



040- Paper Holder

Upholstery Colors



Slate Blue



Wedgewood



Royal Blue



Desert Tan



Burgundy



Country Mist



Warm Gray



Gunmetal



Black