

Family Practice Table with Step Stool

MODEL 8890

- Premium, vinyl upholstery over 2" firm, foam padding (5 cm)
- Full 28.5" width for patient comfort (72 cm)
- Handy, pull-out, welded steel, step stool
- Air-spring adjustable backrest
- Coated, two position adjustable stirrups
- Pull-out leg rest on steel track
- Easy care, Gray, laminate base
- Chip resistant edges
- 1 large side storage compartment with shelf
- Storage compartment is accessible from both sides (no right or left handed tables)
- 2 drawers with all steel slides
- Built with plywood and laminate panels that are California CARB II compliant
- Top sections are easily replaceable in the field
- Adjustable feet for leveling
- Paper dispenser included
- 400 lbs. load capacity under normal use (181 kg)




8890 SPECIFICATIONS


58"/72" extended L x 28.5" W x 31" H



California residents please be advised of the following, pursuant to Proposition 65:

1
 **WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

Features



Laminate Color



Gray

Option



040- Paper Holder

Upholstery Colors



Slate Blue

Wedgewood

Royal Blue

Desert Tan

Burgundy

Country Mist

Warm Gray

Gunmetal

Black