

Panel Leg Series, Space Saver, Exam Table



MODELS: 91100

- All laminate construction
- Sleek, panel leg, space saver design
- Adjustable backrest with 6 settings of approximately: 27°, 36°, 43°, 50°, 55° & 60°
- Pull-out, slate gray, laminate leg rest on steel track extends the table to 72"
- Scratch resistant topcoat finish
- Ships flat for easy handling and extra savings
- Fast and easy assembly
- Seamless rounded corner top
- 2" firm foam padding (5 cm)
- Heavyweight, knit-backed upholstery
- 400 lbs. load capacity under normal use (181.43 kg)

SPECIFICATIONS

91100 58"/72" extended L x 30" W x 31" H




Prop 65 Warning

WARNING: This product can expose you to chemicals, including diisononyl phthalate (DINP), which is known to the State of California to cause cancer. For more information go to www.P65Warnings.ca.gov.

2

Options



034– Paper Dispenser & Cutter Combo

040– Paper Cutter

030– Paper Dispenser

Laminate Colors



Gray



Maple



Dark Cherry

Upholstery Colors



Slate Blue



Wedgewood



Royal Blue



Desert Tan



Burgundy



Country Mist



Warm Gray



Gunmetal



Black