

## Panel Leg Series, Space Saver, Exam Table with Shelf



### MODELS: 91200

- All laminate construction
- Sleek, panel leg, space saver design
- Adjustable backrest with 6 settings of approximately: 27°, 36°, 43°, 50°, 55° & 60°
- Pull-out, slate gray, laminate leg rest on steel track extends the table to 72"
- Scratch resistant topcoat finish
- Ships flat for easy handling and extra savings
- Fast and easy assembly
- Seamless rounded corner top
- 2" firm foam padding (5 cm)
- Heavyweight, knit-backed upholstery
- 400 lbs. load capacity under normal use (181.43 kg)

### SPECIFICATIONS

91100 58"/72" extended L x 30" W x 31" H



### Prop 65 Warning

**WARNING:** This product can expose you to chemicals, including diisononyl phthalate (DINP), which is known to the State of California to cause cancer. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

2

### Options



034- Paper Dispenser & Cutter Combo

040- Paper Cutter

030- Paper Dispenser


### Laminate Colors



Gray



Maple



Dark Cherry

### Upholstery Colors



Slate Blue



Wedgewood



Royal Blue



Desert Tan



Burgundy



Country Mist



Warm Gray



Gunmetal



Black