

Step-Up Table with 4-Doors



MODELS: 9704

- Sturdy, all laminate and steel construction
- Handy, steel stool with slip-resistant top pulls out for use, and neatly stores out of the way when not in use
- Stool can be easily changed on site for use on either side of the table
- Adjustable backrest with 2" for firm foam padding
- 4-door storage accessible from either side of the table
- Non-marking glides for easy positioning
- 400 lbs. weight capacity under normal use

SPECIFICATIONS

9704-27	72 L x 27" W x 31" H
9704-30	72 L x 30" W x 31" H



California residents please be advised of the following, pursuant to Proposition 65:

WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

1

Options



034- Paper Dispenser & Cutter Combo




040- Paper Cutter



030- Paper Dispenser



066- Door and Inside Latch Combo



055- Door Locks



088- Flat Top

Laminate Colors



Dark Cherry



Maple



Gray



Slate Gray

Upholstery Colors



Slate Blue



Wedgewood



Royal Blue



Desert Tan



Burgundy



Country Mist



Warm Gray



Gunmetal



Black