

POWER MAT PLATFORMS



A. MODEL 1440-57-113

A. SERIES 1440 DUAL LIFT POWERMATIC® MAT PLATFORMS




- Power height adjustment from 19" wheelchair access height to 28"H.
- Two direct drive, synchronized lift mechanisms for unsurpassed power.
- Black laminate base with satin aluminum trim and levelers.
- Tear resistant, seamless nylon-reinforced vinyl upholstered top with 2" high-density urethane foam padding for comfort and rounded corners for safety.
- "Hands-free" foot control.
- Hospital grade plug.
- 110V, 60Hz, 2.8 Amps. U.L. Listed


MODEL	WIDTH	LENGTH	HEIGHT
1440-46-113	4'	6'	19"- 28"
1440-46-119	4'	6'	19"- 28"
1440-47-113	4'	7'	19"- 28"
1440-47-119	4'	7'	19"- 28"
1440-57-113	5'	7'	19"- 28"
1440-57-119	5'	7'	19"- 28"
1440-68-113	6'	8'	19"- 28"
1440-68-119	6'	8'	19"- 28"

You must specify upholstery color with order.

UPHOLSTERY - 113 BLUE OR 119 GRAY



SPECIAL ORDER TOP COLORS



B. SERIES 1441 ECONO ELECTRIC HI-LO MAT PLATFORMS



- Power height adjustment from 20" wheelchair height to 29"H.
- Two direct drive, synchronized lift mechanisms for greater power. 600 lbs. patient lift capacity.
- Tear resistant, seamless nylon-reinforced vinyl upholstered top with 2" high-density urethane foam padding for comfort and rounded corners for safety.
- Sturdy black laminate base with leveling legs.
- "Hands-free" foot control.
- Hospital grade plug.
- 110V, 60Hz, 1.3 Amps.

MODEL	WIDTH	LENGTH	HEIGHT
1441-46*	4'	6'	20"-29"
1441-47*	4'	7'	20"-29"
1441-57*	5'	7'	20"-29"

You must specify upholstery color with order.

B. MODEL 1441-47-119

