

# EXAM TABLES



## EXAM TABLES

The medical exam table is essential in any health care facility. Hausted provides high-quality options to meet the needs of hospitals, clinics, ERs and more.

- Five Position Continuous Lock Stirrups, the most positions in the industry
- Pull-Out Leg Rest with Integrated Treatment Tray
  - ♦ Tray is independent of the Legrest and may be accessed easily at any time
- An all-metal construction for the base table provides
- Powder Coat Finish - Foghorn
- Pelvic Tilt



### 4200 Series

**4200-xx** exam table with three pass-through drawers, two storage drawers and stirrups



### 4600 Series

**4600-xx** high/low power exam table with power assisted back, hand control and stirrups. UL listed (60601-1) DC electromechanical drive. Hand pendant controls hi-low, back, flat, return to chair



### Specifications:

|  | 4200 SERIES |        | 4600 SERIES    |                 |
|--|-------------|--------|----------------|-----------------|
|  | inch        | cm     | inch           | cm              |
| Overall Width:                                   | 27.00       | 68.58  | 27.00          | 68.58           |
| Overall Length:<br><i>legrest extended</i>       | 72.00       | 182.88 | 72.00          | 182.88          |
| Overall Height:<br><i>minimum</i>                | 32.00       | 81.28  | 18.50<br>± 0.5 | 46.99<br>± 1.27 |
| Overall Height:<br><i>maximum</i>                | 32.00       | 81.28  | 35.50<br>± 0.5 | 90.17<br>± 1.27 |
| Foot Step Width:                                 | 18.00       | 45.72  | —              | —               |
| Foot Step Length:                                | 17.00       | 43.18  | —              | —               |
| Legrest Width:                                   | 18.00       | 45.72  | 18.00          | 45.72           |
| Legrest Length:                                  | 17.00       | 43.18  | 17.00          | 43.18           |
| Paper Roll Holder                                | 18.00       | 45.72  | 18.00          | 45.72           |
| Widths Available:                                | 19.00       | 48.26  | 19.00          | 48.26           |
|  | 20.00       | 50.80  | 20.00          | 50.80           |
| Drawer Width: <sup>1</sup>                       | 18.00       | 45.72  | —              | —               |
| Drawer Length: <sup>1</sup>                      | 16.50       | 41.91  | —              | —               |
| Drawer Height: <sup>1</sup>                      | 4.00        | 10.16  | —              | —               |
|  |             |        |                |                 |
| Foot Step Capacity:<br><i>EVENLY DISTRIBUTED</i> | lb          | kg     | lb             | kg              |
|  | 800         | 362.87 | —              | —               |
| Weight Capacity:<br><i>EVENLY DISTRIBUTED</i>    | 500         | 226.80 | 600            | 272.16          |
|  | 335         | 151.95 | 395            | 179.17          |
| Shipping Weight:                                 | 335         | 151.95 | 395            | 179.17          |

### Features:

|                     |                               |                               |
|---------------------|-------------------------------|-------------------------------|
| Base:               | pass-through drawer           | dual pedestal                 |
| Stirrups:           | five position continuous lock | five position continuous lock |
| Treatment Tray:     | integrated                    | integrated                    |
| Powder Coat Finish: | foghorn                       | foghorn                       |
| Pull Out Leg Rest:  | standard                      | standard                      |
| Pelvic Tilt:        | standard                      | standard                      |

### Options:

|  | available at additional cost |               |
|--|------------------------------|---------------|
|  | order number                 | order number  |
| None:  | 4200-xx                      | 4600-xx       |
| Drawer Warmer:<br><i>w/ integrated electrical outlet</i> | 4200-DW-xx                   | —             |
| Drawer Warmer*<br><i>includes physician storage</i>      | —                            | 4600-PS-DW-xx |
| Integrated Outlet:                                       | 4200-E0-xx                   | —             |
| Physician Storage*<br><i>without drawer warmer</i>       | —                            | 4600-PS-xx    |
| Assist Arms*   | —                            | 4600-AA-xx    |
| Foot Switch*   | —                            | 4600-FS-xx    |
| Fully Loaded:<br><i>options marked *</i>                 | —                            | 4600-FL-xx    |
| Replacement Top:   | GF420RT-xx                   | GF460RT-xx    |

\* Options Included in the Fully Loaded Package

<sup>1</sup> Interior Dimensions

### Standard Color Options (-xx)

Available in 24 color choices to complement most decors.

|                   |                  |              |                 |                   |                  |
|-------------------|------------------|--------------|-----------------|-------------------|------------------|
| WH<br>Super White | DO<br>Doe        | PM<br>Palm   | SH<br>Shadow    | AR<br>Artichoke   | TE<br>Teal       |
| NO<br>Northwoods  | DB<br>Dusty Blue | PU<br>Plum   | WW<br>Wedgewood | BT<br>British Tan | GR<br>Gray       |
| AL<br>Alabaster   | ST<br>Saddle Tan | PL<br>Pearl  | TL<br>Tea Leaf  | MO<br>Moss        | OG<br>Ocean Gray |
| AQ<br>Aqua Green  | BU<br>Blue       | VI<br>Violet | NA<br>Navy      | WG<br>Warm Gray   | BL<br>Black      |

Swatch Card is available from your representative



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.  
[www.hausted.com](http://www.hausted.com) | [www.grahamfield.com](http://www.grahamfield.com)



© 2018, GF Health Products, Inc. All Rights Reserved.  
Hausted is a registered trademark of GF Health Products, Inc.

Product Hyperlinks Enabled

**graham**  
**field**

GF1800075RevG23