

The image shows a close-up of a white, padded medical examination table. The table is positioned diagonally, with the head of the patient towards the top left. The padding is thick and smooth. The background is a plain, light gray wall. The overall aesthetic is clean and professional.

Ritter®

MEDICAL

Manual Examination Tables

Shop tables & accessories at

RITTER® EXAMINATION TABLES

DESIGNED WITH YOU + YOUR PATIENTS IN MIND

What makes Ritter by Midmark the most purchased brand of examination tables in the industry? It's our willingness to research, study and listen so we can meet your needs—and your budget. The Ritter 204 is the culmination of years of research in design, comfort and function in a manual examination table.



Shop tables & accessories at

2 Manual Examination Tables

Stylish. Comfortable. Easy to clean.

The seamless design of the Ritter 204 upholstered top is as soft and comfortable as it looks. And the seamless upholstery design reduces areas where dirt and bacteria can accumulate so clean-up can be easier, faster and more effective.



REMOVABLE TOP SYSTEM

The one-piece, upholstered top is designed to be easily removed without tools for a more thorough cleaning or a quick upholstery change.



ONE-PIECE SEAMLESS DRAWERS

The Ritter 204 drawers are designed with rounded, seamless edges to help contain spilled fluids and improve infection control.



ENLARGED FOOTSTEP

We increased the size of the footstep to over 40% above the industry standard to help patients feel more secure when accessing and exiting the table.



SPACIOUS SIDE DRAWERS

Three side drawers provide point of care storage for gowns, drapes and other large items, as well as organized storage for smaller supplies and bulk items. Two front drawers provide easy access to items during lower body exams.



STORAGE ACCESS + VOLUME

Pass-through side drawers provide access to supplies from either side of the exam table and give you the flexibility to change your room layout or workflow.



PULLOUT STIRRUPS

Integrated in the exam table, the pullout stirrups provide leg support for patients during lower body examinations.

Accessories + Options that add versatility

Upholstery Options



Flat

Provides an additional working surface for exams and procedures where you need a little more space.



Seamless

Plush, seamless upholstery designed to be more comfortable for the patient and easier to clean.



Soft Touch

Designed to be comfortable for the patient yet durable enough to stand up to the rigors of the clinical environment.

Accessories



Assist Handles

Designed to offer patients a stable platform for accessing the table and a sense of security during ingress, egress and positioning while on the table.



Ritter® 253 LED Exam Light

The 253 is designed with easy access controls, high intensity light and an adjustable focal spot. Available in table or wall mount and mobile options.



Adjustable Arm System

Designed to offer patients a stable platform in accessing the table and a sense of security during positioning. Arms rotate and lock into one of three positions and are easily removable.



Stainless Steel Treatment Pan

Designed for spaces where a stainless steel treatment pan is preferred.



Articulating Knee Crutches

The unique ball-joint design allows for maximum range of adjustment and can help ensure optimum knee and leg support.



Standard Knee Crutches

Swivel to provide leg support for extended procedures where fatigue is a factor or when the patient lacks sufficient leg strength or control.

Ritter® 203

Treatment Table

FEATURES

The economical Ritter 203 is constructed of steel from its frame to its legs. This metal-on-metal design helps give the 203 greater stability and structural integrity. Its flat-top design provides a stable, solid surface for patient manipulations or adjustments during examinations while extra padding and a built-in headrest pillow help keep patients comfortable. And the Mist colored base complements any of the Ritter upholstery colors.

OPTION

Drawers

Rounded drawers with color inserts are available on the right or left side of the table.

EXAMTABLESDIRECT



Ritter® 95 Treatment Table

FEATURES

The Ritter 95 is an economical solution that's built to function in a wide range of applications. It's ideal for new medical practices, industrial first aid areas, school nurses' stations or as an additional table wherever needed. The hardwood construction and H-brace frame design improve strength and stability while the 27" width and 72" length help accommodate a wider range of patients. (Paper roll holder included.)

OPTION

Vinyl Upholstered Shelf

Provides space at the point of care to place supplies or use as an additional worksurface.



Shop tables & accessories at

SPECIFICATIONS

RITTER 204

Height with Upholstery:
Soft Touch Top -
32.13" (81.6 cm)
Seamless Top -
33.13" (84.2 cm)

Patient Weight Capacity:
500 lb (226.8 kg)

Electrical Receptacle:
(115 VAC 5 A Max.)
(230 VAC 2.5 A Max.)
Duplex, hospital grade
electrical receptacle on
left side

Soft Touch Upholstery Top:
30.5" W x 56" L
(77.5 cm x 142.5 cm)

Paper Roll Holder
holds 21" x 4.75" Paper
Roll (53.3 cm x 12 cm)

Length with Footrest
extended: 73.75"
(187.3 cm)

Seamless Upholstery Top:
28.5" W x 59.5" L
(72.4 cm x 151.1 cm)

Paper Roll Holder
holds 21" x 3.5" Paper
Roll (53.3 cm x 8.9 cm)

Length with Footrest
Extended: 76" (193 cm)

UltraFree Upholstery Top:
28.5" W x 59.5" L
(72.4 cm x 151.1 cm)

OPTIONS

Soft Touch, Seamless,
Seamless Ultrafree,
Flat Stitched Upholstery
without pillow,
Flat Ultrafree Upholstery
top without pillow

Pelvic tilt
(1.5° or 7°)

Drawer heater (45 W)*

ACCESSORIES

Knee crutches

Stainless steel
treatment pan

Assist handles

Adjustable arm system

Pass-through drawer stop
(set of 3)

Thick footrest pad kit

ELECTRICAL

115 VAC, 50/60 Hz, 5.0 A
(with receptacle)

115 VAC, 50/60 Hz, 5.5 A
(with drawer heater)

230 VAC, 50/60 Hz, 2.5 A
(with receptacle)

*Not available for
sale in Europe.

RITTER 203

Patient Weight Capacity:
400 lb (182 kg)

Dimensions (top):
27" W x 74" L
(68.6 cm x 188 cm)

Height: 31" (78.7 cm)

Paper Roll Holder:
21" (53.3 cm)

RITTER 95

Patient Weight Capacity:
350 lb (159 kg)

Dimensions (top):
27" W x 72" L
(68.6 cm x 182.9 cm)

Height: 30" (76.2 cm)

Paper Roll Holder:
21" (53.3 cm)

ACCESSORIES



Vinyl Upholstered Shelf

COLOR OPTIONS

STANDARD UPHOLSTERY COLORS

	Mist 813
	Stone 814
	Restful Path 858
	Iron Ore 851
	Lunar Gray 845
	Obsidian 857
	Oasis 843
	Soothing Blue 855
	Dream 847
	Mineral 844
	Healing Waters 848
	Shaded Garden 853
	Citrus 849
	Sandy Retreat 850
	Curative Copper 846
	Cranberry 859
	Dark Linen 856
	Latte 860
	Robust Brown 852
	Deep Earth 854

ULTRAFREE UPHOLSTERY COLORS

	UltraFree Stone 865
	UltraFree Restful Path 868
	UltraFree Obsidian 867
	UltraFree Harbor 861
	UltraFree Wheatgrass 862
	UltraFree Spice 863
	UltraFree Cranberry 869
	UltraFree Dark Linen 866
	UltraFree Latte 870
	UltraFree Branch 864

EXAMTABLESDIRECT



CARB 93129.2 Phase 2 Compliant and TSCA Title VI Compliant

The color examples shown are the best representation of the original material. Actual color may vary slightly. We strongly recommend that you contact Midmark Customer Experience at 1.800.MIDMARK to request a sample before placing your order.

Authorized CE Representative Countries in the EEC should direct all questions, incidents and complaints to the Midmark Authorized CE representative listed below:

CEpartner4U, Esdoornlann 13, 3951 DB Maarn, The Netherlands
Phone: +31.343.442.524 Website: cepartner4u.com

Midmark is an ISO 13485 and ISO 9001 Certified Company. Certain products are not included. See the complete list at: midmark.com/ISO

For more information, contact your Midmark dealer or call:
1.800.MIDMARK
Outside the USA call: 1.937.526.3662 or visit our website: midmark.com

© 2021 Midmark Corporation, Miamisburg, Ohio USA
Products subject to improvement changes without notice. Litho in USA.

007-10261-00 Rev. B1 (9/21)

Shop tables & accessories at

Ritter®