

The image shows a close-up, low-angle view of a white medical examination table. The table is positioned diagonally, with the top right corner of the frame showing the head of the table and the bottom left corner showing the base. The table has a thick, padded white surface. The background is a plain, light-colored wall. The overall lighting is soft and even, highlighting the clean, professional appearance of the equipment.

# Ritter®

MEDICAL

## Manual Examination Tables

# RITTER® EXAMINATION TABLES

## DESIGNED WITH YOU + YOUR PATIENTS IN MIND

What makes Ritter by Midmark the most purchased brand of examination tables in the industry? It's our willingness to research, study and listen so we can meet your needs—and your budget. The Ritter 204 is the culmination of years of research in design, comfort and function in a manual examination table.



# Stylish. Comfortable. Easy to clean.

The seamless design of the Ritter 204 upholstered top is as soft and comfortable as it looks. And the seamless upholstery design reduces areas where dirt and bacteria can accumulate so clean-up can be easier, faster and more effective.



### REMOVABLE TOP SYSTEM

The one-piece, upholstered top is designed to be easily removed without tools for a more thorough cleaning or a quick upholstery change.



### ONE-PIECE SEAMLESS DRAWERS

The Ritter 204 drawers are designed with rounded, seamless edges to help contain spilled fluids and improve infection control.



### ENLARGED FOOTSTEP

We increased the size of the footstep to over 40% above the industry standard to help patients feel more secure when accessing and exiting the table.



### SPACIOUS SIDE DRAWERS

Three side drawers provide point of care storage for gowns, drapes and other large items, as well as organized storage for smaller supplies and bulk items. Two front drawers provide easy access to items during lower body exams.



### STORAGE ACCESS + VOLUME

Pass-through side drawers provide access to supplies from either side of the exam table and give you the flexibility to change your room layout or workflow.



### PULLOUT STIRRUPS

Integrated in the exam table, the pullout stirrups provide leg support for patients during lower body examinations.

# Accessories + Options that add versatility

## Upholstery Options



**Flat**  
Provides an additional working surface for exams and procedures where you need a little more space.



**Seamless**  
Plush, seamless upholstery designed to be more comfortable for the patient and easier to clean.



**Soft Touch**  
Designed to be comfortable for the patient yet durable enough to stand up to the rigors of the clinical environment.

## Accessories



**Assist Handles**  
Designed to offer patients a stable platform for accessing the table and a sense of security during ingress, egress and positioning while on the table.



**Ritter® 253 LED Exam Light**  
The 253 is designed with easy access controls, high intensity light and an adjustable focal spot. Available in table or wall mount and mobile options.



**Adjustable Arm System**  
Designed to offer patients a stable platform in accessing the table and a sense of security during positioning. Arms rotate and lock into one of three positions and are easily removable.



**Stainless Steel Treatment Pan**  
Designed for spaces where a stainless steel treatment pan is preferred.



**Articulating Knee Crutches**  
The unique ball-joint design allows for maximum range of adjustment and can help ensure optimum knee and leg support.



**Standard Knee Crutches**  
Swivel to provide leg support for extended procedures where fatigue is a factor or when the patient lacks sufficient leg strength or control.

# Ritter® 203 Treatment Table

## FEATURES

The economical Ritter 203 is constructed of steel from its frame to its legs. This metal-on-metal design helps give the 203 greater stability and structural integrity. Its flat-top design provides a stable, solid surface for patient manipulations or adjustments during examinations while extra padding and a built-in headrest pillow help keep patients comfortable. And the Mist colored base complements any of the Ritter upholstery colors.

## OPTION

### Drawers

Rounded drawers with color inserts are available on the right or left side of the table.



# Ritter® 95 Treatment Table

## FEATURES

The Ritter 95 is an economical solution that's built to function in a wide range of applications. It's ideal for new medical practices, industrial first aid areas, school nurses' stations or as an additional table wherever needed. The hardwood construction and H-brace frame design improve strength and stability while the 27" width and 72" length help accommodate a wider range of patients. (Paper roll holder included.)

## OPTION

### Vinyl Upholstered Shelf

Provides space at the point of care to place supplies or use as an additional workspace.



Shop tables & accessories at



# SPECIFICATIONS

## RITTER 204

Height with Upholstery:  
Soft Touch Top -  
32.13" (81.6 cm)  
Seamless Top -  
33.13" (84.2 cm)

Patient Weight Capacity:  
500 lb (226.8 kg)

Electrical Receptacle:  
(115 VAC 5 A Max.)  
(230 VAC 2.5 A Max.)  
Duplex, hospital grade  
electrical receptacle on  
left side

Soft Touch Upholstery Top:  
30.5" W x 56" L  
(77.5 cm x 142.5 cm)

Paper Roll Holder  
holds 21" x 4.75" Paper  
Roll (53.3 cm x 12 cm)

Length with Footrest  
extended: 73.75"  
(187.3 cm)

Seamless Upholstery Top:  
28.5" W x 59.5" L  
(72.4 cm x 151.1 cm)

Paper Roll Holder  
holds 21" x 3.5" Paper  
Roll (53.3 cm x 8.9 cm)

Length with Footrest  
Extended: 76" (193 cm)

UltraFree Upholstery Top:  
28.5" W x 59.5" L  
(72.4 cm x 151.1 cm)

## OPTIONS

Soft Touch, Seamless,  
Seamless Ultrafree,  
Flat Stitched Upholstery  
without pillow,  
Flat Ultrafree Upholstery  
top without pillow

Pelvic tilt  
(1.5° or 7°)

Drawer heater (45 W)\*

## ACCESSORIES

Knee crutches

Stainless steel  
treatment pan

Assist handles

Adjustable arm system

Pass-through drawer stop  
(set of 3)

Thick footrest pad kit

## ELECTRICAL

115 VAC, 50/60 Hz, 5.0 A  
(with receptacle)

115 VAC, 50/60 Hz, 5.5 A  
(with drawer heater)

230 VAC, 50/60 Hz, 2.5 A  
(with receptacle)

\*Not available for  
sale in Europe.

## RITTER 203

Patient Weight Capacity:  
400 lb (182 kg)

Dimensions (top):  
27" W x 74" L  
(68.6 cm x 188 cm)

Height: 31" (78.7 cm)

Paper Roll Holder:  
21" (53.3 cm)

## RITTER 95

Patient Weight Capacity:  
350 lb (159 kg)

Dimensions (top):  
27" W x 72" L  
(68.6 cm x 182.9 cm)

Height: 30" (76.2 cm)





















Paper Roll Holder:  
21" (53.3 cm)

## ACCESSORIES

Vinyl Upholstered Shelf

# COLOR OPTIONS

## STANDARD UPHOLSTERY COLORS

	Mist 813
	Stone 814
	Restful Path 858
	Iron Ore 851
	Lunar Gray 845
	Obsidian 857
	Oasis 843
	Soothing Blue 855
	Dream 847
	Mineral 844
	Healing Waters 848
	Shaded Garden 853
	Citrus 849
	Sandy Retreat 850
	Curative Copper 846
	Cranberry 859
	Dark Linen 856
	Latte 860
	Robust Brown 852
	Deep Earth 854

## ULTRAFREE UPHOLSTERY COLORS

	UltraFree Stone 865
	UltraFree Restful Path 868
	UltraFree Obsidian 867
	UltraFree Harbor 861
	UltraFree Wheatgrass 862
	UltraFree Spice 863
	UltraFree Cranberry 869
	UltraFree Dark Linen 866
	UltraFree Latte 870
	UltraFree Branch 864

# EXAMTABLES DIRECT



CARB 93129.2 Phase 2 Compliant and TSCA Title VI Compliant

The color examples shown are the best representation of the original material. Actual color may vary slightly. We strongly recommend that you contact Midmark Customer Experience at 1.800.MIDMARK to request a sample before placing your order.

Authorized CE Representative Countries in the EEC should direct all questions, incidents and complaints to the Midmark Authorized CE representative listed below:

CEpartner4U, Esdoornlann 13, 3951 DB Maarn, The Netherlands  
Phone: +31.343.442.524 Website: [cepartner4u.com](http://cepartner4u.com)

Midmark is an ISO 13485 and ISO 9001 Certified Company. Certain products are not included. See the complete list at: [midmark.com/ISO](http://midmark.com/ISO)

For more information, contact your Midmark dealer or call:  
1.800.MIDMARK  
Outside the USA call: 1.937.526.3662 or visit our website: [midmark.com](http://midmark.com)

© 2021 Midmark Corporation, Miamisburg, Ohio USA  
Products subject to improvement changes without notice. Litho in USA.

007-10261-00 Rev. B1 (9/21)

# Ritter®

Shop tables & accessories at