



The Sonesta 6202 is the ultimate in positioning for your practice. Popular accessories, such as fold up foot rests and a seat cushion without cut out are included standard with model 6202. Additionally the locking mechanism for one touch movement, offers best-in-class comfort for the patient and the best available ergonomics and ease-of-use for the physician.

## Features

- 5 functions, maneuvered with 3 motors (height, back and seat position) standing position.
- Adjusted by back motor, maneuvered using the hand control.
- Rotation of seat/back with locked position.
- Adjusted by back motor, maneuvered with hand control.



# Specifications

## General Specifications

CE-marked in accordance with guidelines in Medical Device Directive 93/42EEC.

Length without foot rest: 156cm / 63".

Length without foot rest: 196cm / 78".

Width (total): 99cm /39".

Height with seat cushion (min.): 50cm /20".

Height with seat cushion (max.): 130cm /51".

Patient weight: max.250kg/551lbs.

Wheel diameter: 12,5cm/4".

Net weight: Approx 150 kg / 330 lbs.

Power consumption: max 5.8A (120V~), 3.0A (230V~).

Input voltage: 120/230V ~ 50/60 Hz.

Output voltage: 24v DC.

Satisfies the requirements for EMC according to IEC 60601-1-2.