

# Fusion**ONE+**<sup>TM</sup>

Power Hi-Lo, Power Back Exam Table with Foot Control, Drain Pan, Stirrups, Drawer Warmer, Pelvic Tilt & Receptacle

Model 3503



## For examinations

Fusing the accessibility of a power table and the affordability of a fixed-height box table, the FusionONE+ Power Exam Table offers the features healthcare providers need most.

Featuring a power back, this hi-lo exam table can accommodate multiple specialties. It enables simple ADA access and ease of examination for patients of all ages and abilities. The power backrest makes it easy to lie patients back and sit them up safely and comfortably.

**PROTECTS** > patients from infection > facilities from damage > staff from injury

## Key features at a glance

- > Power hi-lo and power backrest
- > Two-function foot control
- > 18"-21" (45.7 cm-53.3 cm) paper roll holder
- > 500 lb. (226.8 kg) patient weight capacity
- > ADA Accessible 19" (48.3 cm)
- > Pelvic tilt (7°)
- > Drain pan: plastic 1.5" deep (3.8 cm)
- > Large, reversible side drawer and large front drawer
- > Removable, seamless upholstered top with PreFixx<sup>®</sup> Protective Finish
- > Durable, galvanized steel construction
- > Electrical receptacle
- > Drawer warmer
- > Heavy-gauge steel pull-out leg extension
- > Foldaway stirrups with four lateral positions and adjustable length
- > Easy-clean powder coat enamel surfaces
- > 3-year limited warranty

## Specifications

- › Height: min. 19" (48.3 cm) at seat max. 32" (81.3 cm)
- › Width: 28" (71.1 cm)
- › Length: 58"-76" (147.3 cm-193 cm) with leg rest extended
- › Leg rest: 18" W x 18" D (45.7 cm W x 45.7 cm D)
- › Storage Drawer (Front): 18" W x 12" D x 3" H (45.7 cm x 30.5 cm x 7.6 cm)
- › Storage Drawer (Side): 17.5" W x 21" D x 5" H (44.5 cm x 53.3 cm x 12.7 cm)
- › Single-column, electro-magnetic actuator
- › Electrical: 120 VAC 10.5 AMP / 220 VAC 5.8 AMP, 50-60Hz
- › Weight capacity: 500 lbs. (226.8kg)
- › Net weight: 350 lbs. (159kg)
- › Shipping weight: 390 lbs. (177kg)
- › Shipping dimensions: 30" W x 64" L x 31" H (76.2 cm x 162.6 cm x 78.7 cm)

## Accessories

- › 251 – Bierhoff Knee Crutches (pair)
- › 582 – Hand Control
- › 585 – Armrests
- › 590 – Welch Allyn Light Bracket capabilities for #48955 (patient left)
- › 591 – Welch Allyn Light Bracket capabilities for #48955 (patient right)

## Colors



**ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA**

©2023 United Metal Fabricators Inc dba UMF Medical 04.23 All Rights Reserved. PreFixx® is a registered trademark of Omnova Solutions. FusionONE™ is a trademark of United Metal Fabricators.