

Bariatric & Parent with Child Phlebotomy Chair

Model 8690

UMF Medical's 8690 Bariatric Blood Draw, "Mommy-and-Me" Chair offers a flip-up, height-adjustable arm board & front board. Heavy-duty steel frame, with a 750lb weight capacity. Available in all 16 upholstery colors with PreFixx® Protective Finish.



Colors



ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA



Key features at a glance

- › Ultra comfort seamless upholstery
- › Accommodates bariatric patients and parent with child
- › Padded, flip-up, height adjustable arm board
- › Heavy duty, all-welded steel frame
- › 750lb (340.2kg) patient weight capacity
- › Adjustable leg levelers
- › 3-year limited warranty

Specifications

Width: 41" (104cm)
Depth: 35" (88.9cm)
Height: 41.75" (106cm)
Net Weight: 78lb (35.4kg)
Shipping Weight: 138lb (62.6kg)
Shipping Dimensions: L44xW34xH50
(112cm x 86.4cm x 127cm)

Seat dimension:
Height: 22" (55.9cm)
Width: 26.625" (68cm)
Between arm boards:
30" (76cm)

Accessories

120 – Side tray
121 – Storage Box